Performance of a novel NGS-based Bridge Capture™ technology for mutation profiling in liquid biopsies derived from metastatic colorectal cancer patients

genomill UNIVERSITY OF TURKU

Anttoni Korkiakoski^{1,2}, Simona Adamusová^{1,2}, Aparna Ganesan^{3,4}, Anna Musku¹, Tuula Rantasalo¹, Nea Laine¹, Emma Andersson⁵, Emerik Österlund^{6,7}, Ali Ovissi⁸, Päivi Halonen⁹, Tatu Hirvonen¹, Jorma Kim¹, Jukka Laine^{1,10}, Antti Silvoniemi^{10,11}, Heikki Minn¹⁰, Juuso Blomster¹, Anna-Kaisa Anttonen⁵, Soili Kytölä⁵, Pia Osterlund^{9,12,13}, Juha-Pekka Pursiheimo¹, Pirjo Nummela^{3,4}, Manu Tamminen^{1,2} & Ari Ristimäki^{3,4}

Background

Liquid biopsy is a minimally invasive alternative to tissue biopsy and imaging modalities for cancer diagnostics. Circulating tumor DNA (ctDNA), tumor-derived DNA in the bloodstream, can be used to monitor treatment response, detect minimal residual disease (MRD), and identify resistance-associated mutations. This is the first study to evaluate the performance of a novel NGS-based Bridge Capture™ technology for mutation profiling and MRD detection in ctDNA from patients with metastatic colorectal cancer (mCRC).

Study Design

The performance of Bridge Capture™ technology was compared to droplet digital PCR (ddPCR), Ion AmpliSeq Cancer Hotspot Panel v2 (CHPv2), and RT-qPCR based Idylla ctKRAS Mutation Assay. Additionally, the analytical performance and scalability of the technology were assessed using synthetic DNA targets.

Serial plasma samples from ten patients with KRAS mutation in the primary tumor tissue (n = 80)




Figure 2. Evaluation of Bridge Capture™ using (a) ctDNA extracted from plasma samples of mRCR patients and (b) synthetic DNA targets. ddPCR: droplet digital PCR; NGS: next generation sequencing; RCA: rolling circle amplification; VAF: variant allele frequency

Bridge Capture™




Figure 1. Bridge Capture™ single-tube workflow is simple and rapid. It allows first-step sample indexing, and after only few steps libraries for new generation sequencing platforms can be produced. Red circle depicts a mutation, purple, pink, and orange squares depict a sample index, and pink circle is a F

Results

Bridge Capture™ showed substantial agreement and perfect agreement with ddPCR (Cohen's Kappa 0.70), and Idylla (Cohen's Kappa 0.79), and with Ion AmpliSeq CHPv2, respectively (Table 1). Bridge Capture™ displayed very strong correlation between both ddPCR (ρ = 0.86) (Figure 3a) and Ion AmpliSeq ($\rho = 0.74$) (Figure 3b). Bridge Capture™ expanded 851-probe panel showed very strong correlation until 0.03% (R = 0.99)(Figure 3c). Figure 4 displays serial samples of four mCRC patients with KRAS mutations detected by ddPCR and Bridge Capture™, including eight additional CRC driver mutations in other genes mirroring the KRAS mutations.

Table 1. Concordance between Bridge Capture[™] and ddPCR, Idylla, and Ion AmpliSeq CHPv2 in detecting *KRAS* mutations from ctDNA of mCRC patients. Cohen's Kappa measures the agreement between two raters. The range between 0.61–0.8 is regarded as substantial agreement.

		Bridge Capture™		Cohen's
		+	_	Kappa
ddPCR	+	35	4	0.70
(n = 80)	_	8	33	
ldylla	+	31	4	0.79
(n = 58)*	_	2	21	
lon AmpliSeq	+	7	0	1.0
	_	0	3	

^{*7} samples failed (1 positive and 6 negative for Bridge Capture™).




Figure 3. Correlation of ctDNA VAF values from mCRC patients between Bridge Capture™ and other methods. (a) Comparison with ddPCR for KRAS mutations across 80 patient samples (ρ = 0.86). (b) Comparison with Ion AmpliSeq CHPv2 across 10 patient samples (ρ = 0.74). (c) Dilutions from 0.0313–1% displayed strong linear correlation between the observed and expected total specific VAF values (Pearson R = 0.99). Violet circles are mean values of 5 replicates and SD is displayed as error bars. Individual linear regressions of each probe are displayed as violet lines.




Figure 4. Serial samples of four mCRC patients with KRAS (circle) mutations detected by Bridge Capture™ (violet) and ddPCR (red). Bridge Capture™ reported additional oncogenic mutations (gray) in TP53 (square), APC (diamond), and PIK3CA genes (triangle).

Conclusions

technology compatible with all new generation sequencing platforms. The technology's unique features enable first-step sample indexing which simplifies even further the already fewstep single-tube sample processing enabling increased cost benefits.

Bridge Capture™ is simple, rapid and Bridge Capture™ improved mutation profiling sensitive targeted sequencing and MRD detection in mCRC ctDNA samples by reporting the broader mutational landscape without sacrificing sensitivity. It displayed substantial agreement with other technologies and the observed scalability enables even more comprehensive view of the mutational landscape without the need for extremely deep sequencing.

¹Genomill Health Inc, Turku, Finland; ²University of Turku, Department of Biology, Turku, Finland; ³University of Helsinki, Applied Tumor Genomics Research Program, Helsinki, Finland; ⁴Helsinki University Hospital and University of Helsinki, Department of Pathology, Helsinki, Finland; ⁵Helsinki University Hospital and University of Helsinki, Department of Genetics, Helsinki, Finland; ⁶Helsinki University Hospital and University of Helsinki, Department of Surgery, Helsinki, Finland; ⁷Uppsala University, Department of Immunology, Genetics and Pathology, Uppsala, Sweden; ⁸Helsinki University Hospital and University of Helsinki, Department of Radiology, Helsinki, Finland; ⁹Helsinki University Hospital and University of Helsinki, Department of Oncology, Helsinki, Finland; 10Turku University Hospital, Department of Oncology, Turku, Finland; 11Turku University Hospital and University of Turku, Department of Otorhinolaryngology, Turku, Finland; 12Tampere University Hospital and University of Tampere, Department of Oncology, Tampere, Finland; 13Karolinska University Hospital and Karolinska Institute, Department of Oncology/Pathology, Stockholm, Sweden