

GENO1®: A NOVEL CIRCULAR LIBRARY PREP FOR RAPID AND UNIVERSAL NGS ANALYSIS

Simona Adamusová^{1,2}, Anttoni Korkiakoski^{1,2}, Nea Laine¹, Aparna Ganesan^{3,4}, Tatu Hirvonen¹, Anna Musku¹, Tuula Rantasalo¹, Jorma Kim¹, Juuso Blomster^{1,2}, Jukka Laine^{1,2,5}, Ian McLaughlin⁶, Pirjo Nummela^{3,4}, Ari Ristimäki^{3,4}, Manu Tamminen^{1,2}, Juha-Pekka Pursiheimo¹

¹Genomill Health Inc., Turku, Finland; ²University of Turku, Turku, Finland; ³University of Helsinki, Applied Tumor Genomics Research Program, Helsinki, Finland; ⁴Helsinki University Hospital and University of Helsinki, Department of Pathology, Helsinki, Finland; ⁵Department of Pathology, TYKS Laboratories, Turku University Hospital, Turku, Finland; ⁶Pacific Biosciences, Menlo Park, CA, USA

LIBRARY PREP – BOTTLENECK OF LIQUID BIOPSIES

Traditional library prep methods are **slow, labor-intensive and complex**, increasing operational demands and negatively impacting decentralised laboratories.

Unnecessary complexity introduces errors and reduces data quality, thereby limiting clinical utility.

GENO1® REMOVES THE BOTTLENECK

RAPID & SIMPLE

Streamlined library prep from the very first step, with **minimal hands-on time and fast turnaround**. Eliminates the need for specialized equipment, supporting use in diverse settings.

EXTREMELY SENSITIVE

Detects rare mutations, enabling early cancer detection and minimal residual disease monitoring.

COST-EFFICIENT

Simplified workflow and integration of early sample indexing lower overall costs.

PLATFORM AGNOSTIC


The PCR-free workflow generates multiple library formats, allowing seamless integration across a wide range of NGS platforms.

FIND FULL DETAIL




IN RECENT STUDIES

1 DNA CIRCULARIZATION WITH EARLY SAMPLE INDEXING



Twenty-five different Loops competed for their incorporation into Converted Circular DNA. The Loops performed comparably with **no preferential enrichment**, demonstrating feasibility of early sample indexing.




4 CLINICAL VALIDATION

Serial plasma samples (n = 80) from ten colorectal cancer patients with known KRAS variant were analyzed. **Geno1®** showed **substantial and perfect agreement** with ddPCR and Idylla, and with Ion AmpliSeq CHPv2, respectively (Cohen's Kappa).




	GENO1®		Cohen's Kappa
	+	-	
ddPCR (n = 80)	35	4	0.70
	8	33	
Idylla (n = 58)	31	4	0.79
	2	21	
Ion AmpliSeq (n = 10)	7	0	1.0
	0	3	

Matched KRAS variant detected by **Geno1®** (violet) and by reference method **ddPCR** (red). **Geno1®** revealed also **previously undetected variants** (grey), analysed with a 282-probe panel.



2 DIRECT SEQUENCING OF CIRCULAR LIBRARIES

Direct circular sequencing on PacBio Onso™ yielded **high Q scores** and **uniform read composition** across ten probes for both forward and reverse reads.



3 SHORT-READ NANOPORE SEQUENCING

Concatemers composed of short-read repeats form Nanopore-compatible libraries. Alignment of repeats into consensus sequence improves accuracy with repeat count, reaching unprecedented Q scores (up to 1000).

